

# An observational cohort study of Multiplex8+: a spatially informed assay that uses multiplexed RNA-FISH guided laser capture microdissection followed by total RNA-sequencing

Evan D. Paul<sup>1,2,\*</sup>, Barbora Huraiová<sup>1,2</sup>, Natália Matyášová<sup>1,2,3</sup>, Daniela Gábrišová<sup>1,2</sup>, Soňa Gubová<sup>1,2</sup>, Tomáš Ondříš<sup>1,2</sup>, Michal Gala<sup>1,2</sup>, Liliane Barroso<sup>1,2</sup>, Helena Ignačáková<sup>1,2</sup>, Natália Valková<sup>1,2</sup>, Fresia Pareja<sup>4</sup>, Jakob N. Kather<sup>5,7</sup>, Pavol Čekan<sup>1,2</sup>

<sup>1</sup> MultiplexDX, s.r.o., Comenius University Science Park, Bratislava, Slovakia.; <sup>2</sup> MultiplexDX, Inc., Rockville, MD, USA.; <sup>3</sup> Institute of Clinical Biochemistry and Diagnostics, University Hospital, Faculty of Medicine in Hradec Králové, Charles University, Hradec Králové, Czech Republic; <sup>4</sup> Department of Pathology and Laboratory Medicine, Memorial Sloan Kettering Cancer Center, New York, NY, USA. <sup>5</sup> Else Kroener Fresenius Center for Digital Health, Technical University Dresden, Dresden, Germany.; <sup>6</sup> Department of Medicine I, University Hospital Dresden, Dresden, Germany.; <sup>7</sup> Medical Oncology, National Center for Tumor Diseases (NCT), University Hospital Heidelberg, Heidelberg, Germany. \* Presenting author.


**P1-05-23**

## Background

- On a retrospective cohort of 1,082 FFPE breast tissues, we previously developed and validated mFISHseq, a novel, spatially informed tool that integrates multiplexed RNA fluorescent in situ hybridization (FISH) of the four main breast cancer biomarkers (*ESR1/PGR/ERBB2/MKI67*), which are used to guide laser capture microdissection (LCM) followed by RNA-sequencing<sup>1-4</sup>.
- Here, we applied a Research Use Only (RUO) version of this test called Multiplex8+ to an observational cohort of 53 patients, where we assessed intrinsic molecular and TNBC subtypes, prognostic risk, and the expression of 40 genes and 28 gene signatures related to cancer hallmark pathways and treatment response.
- Objectives:** We compared Multiplex8+ to immunohistochemistry and associated our list of genes and gene signatures with response to a variety of treatments, including chemotherapies and novel targeted therapies like CDK4/6 inhibitors and antibody drug conjugates (ADCs).

## Methods

### The Multiplex8+ assay<sup>1</sup>



- ACD RNAScope Multiplex Fluorescent V2 Assay
- Estrogen (*ESR1*)
- Progesterone (*PGR*)
- Her2 (*ERBB2*)
- Ki67 (*MKI67*)
- H&E + Biomarker-guided capture of ROIs
- Takara SMARTer Stranded Total RNA-Seq Kit v3 - Pico Input Mammalian
- NovaSeq 6000 -100M reads/sample (2 x 100 PE)
- 4 key biomarkers
- Molecular/TNBC subtype
- Prognostic risk
- Genes and gene signatures
- Recommendations

### Genes / gene signatures

Drug type / pathway	Genes	Gene signatures
Prognosis	1	3
Proliferation	2	2
Luminal	2	4
Her2	3	3
Chemotherapy	9	8
Immunotherapy	3	5
CDK4/6 inhibitors	5	0
PI3K inhibitors	0	1
DNA damage and repair	2	1
Angiogenesis / Hypoxia	1	1
Antibody drug-conjugates	20	0
Total	40	28


### Overview of RUO cohort

Characteristic	HR+/HER2-	HR+/HER2+	HR-/HER2+	HR-/HER2-	Total
No. of patients	10	9	3	30	53*
Average age (years)	51	49	54	45	47
<b>Tumor size (T)</b>					
pT1	3	3	2	12	20 (38%)
pT2	4	1	0	5	10 (19%)
pT3 + pT4	1	2	0	3	6 (11%)
Unavailable	2	3	1	10	17 (32%)
<b>Node status (N)</b>					
Negative	3	1	1	14	19 (36%)
Positive	4	4	0	5	13 (25%)
Unavailable	3	4	2	11	21 (40%)
<b>Grade (G)</b>					
G1	1	1	0	0	2 (4%)
G2	2	5	0	8	15 (28%)
G3	5	1	3	18	27 (51%)
Unavailable	2	2	0	4	9 (17%)


\* One patient did not have clinicopathological information.

References: 1. Paul, E.D. et al. medRxiv 2023.12.05.23299341, doi: <https://doi.org/10.1101/2023.12.05.23299341> (2023). 2. Paul, E.D. et al. Annals of Oncology (2024) 9 (suppl\_4), 1-34. 10.1016/esoop/esoop103010. 3. Loeffler, C.M.L. et al. J Clin Oncol 42, 2024 (suppl 16; abstr 3069). 4. Pareja F. et al. Annals of Oncology (2024) 35 (suppl\_2): S309-S348. 10.1016/annonc/annonc1577.


## Results – IHC concordance and predictive efficacy



**2** Multiplex8+ has high predictive accuracy, especially for targeted therapies, where 91% (10/11) of the predictions were correct as supported by patient follow-up data regarding clinical response. By combining ADC-processing biomarkers, Multiplex8+ may stratify patients into ADC-treatment responsive subgroups.



## Results – Patient vignettes



## Summary and Conclusions

- We demonstrated the feasibility of implementing Multiplex8+ in a real-world setting and illustrated its clinical potential in predicting response to a variety of treatments.
- Multiplex8+ showed high concordance with IHC (91%), mirroring the 93% accuracy observed in our large-scale retrospective cohort of 1,082 breast tissues<sup>1</sup>.
- By interrogating the expression of targets related to ADC antigens, cytotoxic payloads, endocytosis, lysosome function, and resistance, we illustrate a hypothetical framework that could be used in a future prospective validation for identifying patients responsive to ADCs.

## Funding and contact

Funding: MultiplexDX and the European Union's Horizon 2020 research and innovation programme under an EIC Accelerator grant (agreement No 946693).

Scan here to learn more and download the poster



Copies of this poster obtained through Quick Response (QR) Code are for personal use only and may not be reproduced without permission from SABCS® and the author of this poster.

Contacts: paul@multiplexdx.com  
pavol@multiplexdx.com